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If you want other patient you must give a different ID.Do you really want to add this patient in your folder?

YES **NO**

Patient ID	123450
Nat. Healthcare ID	
First Name	Andreas
Last Name	Savva
Gender	
Date of Birth	Aug. 13, 1964
Address	Grigori Afentiou
Street No.	42
Post code	8512
Province	Kamaries
City	Larnaca
Country	Cyprus
Telephone	99999999
Email	

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Addition of a New Patient

Patient's Information

Patient ID*		Nat. Healthcare ID	
First Name*		Last Name*	
Birth Date*	/ /	Gender*	Male
Address		Street No.	
Province		Post code	
City		Country	
Telephone*		Email	

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Doctor's Information	
Doctor ID	123456
National Doctor ID	
First Name	Stavros
Last Name	Stamatopoulos
Gender	Male
Date of Birth	Nov. 6, 1972
Address	
Street No.	6
Post code	
Province	
City	Larnaca
Country	Cyprus
Telephone	99999999
Email	stamatop.stavros@gmail.com

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Doctor's Information

Doctor ID*	123456	National Doctor ID	
First Name*	Stavros	Last Name*	Stamatopoulos
Birth Date*	8/11/1972	Gender*	Male
Address	Axeaspolis	Street No.	56
Province		Post code	
City	Larnaca	Country	Cyprus
Telephone*	9999999	Email	stamatop.stavros@gmail.com

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
A screenshot of a web application interface for editing a patient's profile. The title bar at the top says "Edit the patient's profile". Below it, there's a section titled "Patient's Information". The form contains two columns of input fields. The left column includes: Patient ID* (123450), First Name* (Andreas), Birth Date* (13/8/1964), Address (Grigori Alevtiov), Province (Kamaries), City (Larnaca), and Telephone* (99999999). The right column includes: Nat. Healthcare ID, Last Name* (Savva), Gender* (Female), Street No. (42), Post code (8512), Country (Cyprus), and Email. A stethoscope is visible in the background on the right side of the form.

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addressData:

- Address_id (Πρωτεύον κλειδί)
- Street
- City
- Post_code
- Province
- Country
- Street_no

doctorData:

- Di (Πρωτεύον κλειδί)
- Doctor_id
- National_healthcare_doctor_id
- Given_name
- Family_name
- Date_of_birth
- Gender_coding
- Address
- Telephone_no
- Email
- Pass_word

patientData:

- Id (Πρωτεύον κλειδί)
- Patient_id
- National_healthcare_patient_id
- Given_name
- Family_name
- Date_of_birth
- Gender_coding
- Address
- Telephone_no
- Email

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
Operation Manual

Sign up

For a doctor to use the system, a registered account is required.
Figure 1 – Admin log in page

A screenshot of the admin login page. At the top, a blue banner contains the text "Welcome Admin to the Online System for Ultrasound Carotid Analysis!" and "LOGIN to continue" in red. Below the banner is a "Sign In" form with two input fields for "Email" and "Password", a blue "Login" button, and a link "Forgot password?". The background of the page features a light blue gradient with a stethoscope graphic on the right side.

An account can be created only by an admin, to avoid any problems and to ensure the confidentiality of the medical data, in case someone tries to use the system without being a doctor.

To sign up a doctor, the admin must sign to the link http://173.249.18.105:8000/admin_login/ and fill the form (Figure 2) with the new doctor's data. The fields marked with a start (*) must be filled and in case admin tries to submit the form while one of them is left blank, then a message alert (Figure 3) pops up to inform the admin about filling them. Press the second button at the left-up part of the page to save doctor's information.


A screenshot of the "Addition of a New Doctor" form. The form is titled "Doctor's Information" and is divided into two columns. The left column contains fields for "Doctor ID*", "First Name*", "Birth Date*" (with a date picker), "Address", "Province", "City", and "Telephone*". The right column contains fields for "Nat. Health. Dr ID", "Last Name*", "Gender*" (with a dropdown menu), "Street No.", "Post code", "Country", and "Email". The background of the page features a light blue gradient with a stethoscope graphic on the right side.

Figure 2 – Form to add a new doctor to the system. Fields with the symbol * are necessary to be filled. By pressing the first button at the left up part of the page admin goes to previous page. By pressing the second button, doctor's information are saved in the database. By pressing the third button, page refreshes, in case a mistake is made.


Figure 3 – Message alert in case necessary filled fields(*) are left blank.

Doctor's Sign in to the System

To use the system you must first log in. Go to the <http://173.249.18.105:8000/login>, where the log in page (Figure 4) will show up.

Figure 4 – Log in page

Use your email account and your secret password to log in and you would be transferred to <http://173.249.18.105:8000/menu> page. If the email is not valid a message appears to inform you (Figure 5).


Figure 5 – Message about invalid information during log in.

In case you forget your password you can choose “Forgot Password” (Figure 4, bottom-center part of the page). You would be transferred at http://173.249.18.105:8000/forgot_pass.

New Password

You can get a new password at http://173.249.18.105:8000/forgot_pass (Figure 6).



Figure 6 – Forgot your password page.

In the input box, write your email that is registered to the system. If the email is invalid, then a message alerts you (Figure 7).

Figure 7 – Message in case of invalid email

If the email is valid, a new, random password is sent to your email (Figure 8) from the account epl400ade2018@gmail.com and you are transferred at the log in page (Figure 3). After, you can use your email and your new password to log in.


Figure 8 – Email with a new password.

Main Menu

After Sign in user goes to Main Menu (Figure), at <http://173.249.18.105:8000/menu>.


Figure 9 – Main Menu

Page gives the following options:

6. Patients
7. Archive
8. Manual
9. Settings
10. Exit

Patients

By choosing “Patients”, you can choose between option “Add Patient” or “Browse Patient” (Figure 10)



Figure 10 – Options “Add Patient” and “Browse Patient”

Option “Add Patient” leads to page http://173.249.18.105:8000/add_new_patient/. An input dialog asks for new patient’s ID (Figure 11).


Enter the ID of the new patient.

OK Cancel

Figure 11 – Input dialog for patient’s ID.

In case patient exists in database from another doctor, then a message alert shows up with all patient’s personal information to verify (Figure 12).

The ID already exists. The details of the corresponding patient is shown below.


If you want other patient you must give a different ID. Do you really want to add this patient in your folder?

YES **NO**

Patient ID	123450
Nat. Healthcare ID	
First Name	Andreas
Last Name	Savva
Gender	
Date of Birth	Aug. 13, 1964
Address	Grigori Afventiou
Street No.	42
Post code	8512
Province	Kamaries
City	Larnaca
Country	Cyprus
Telephone	99999999
Email	

Figure 12 – Message to inform user that patient already exists.

If patient is new to the system, a form is shown to fill up with patient’s personal information (Figure 13).

Navigation icons: back, forward, search, etc.

Addition of a New Patient

Patient's Information

Patient ID*		Nat. Healthcare ID	
First Name*		Last Name*	
Birth Date*	/ /	Gender*	
Address		Street No.	
Province		Post code	
City		Country	
Telephone*		Email	

Figure 13 – Patient registration form. Fields with * must be filled.

Fields marked with star (*) must be field. If any of them is left black and user tries to submit the form, a message alert warns him to fill the rest of the fields.

With option “Browse patient”, user can browse a patient that he already added to the system. When choosing this button, a window (Figure 14) with all available patients appears.

Figure 14 – Doctor’s available patients.

After choosing a patient, user goes to the Operations page (Figure 24).

Archives

The button “Archives” gives you the options “IMT data” and “Plaque data”. In the first option, you can see the Intima Media Thickness measurements report while in the second option you can see the Plaque measurements report.

IMT measurements report

If user chooses the option “IMT Report”, user can choose the patient he wants from a list of all patients he has (Figure 14). To continue, user must also choose from a list the date and time of all previous reports of the same patient (Figure 15).


Figure 15 – List of previous measurements.

On the screen a report with all previous measurements taken at the Define Image Scale/Measurements operation is shown on screen (Figure 16).


Figure 16 –IMT Report

By selecting Print Results button user print the results from his personal printer.

By selecting Save Results user can save results in his personal computer with the name “IMT_Results.xlsx”.

By selecting Close user goes back to Main Menu (Figure 9).

Plaque measurements report

A pop-up window showing all the available patients appears. You must select the patient you want to see the measurements of plaque, then another pop-up window showing the available datetimes when there are metrics recorded on the database, is on the screen.

When you make a choice, you can see all the plaque measurements, for the corresponding patient, on the date and time you chose (Figure 17).


Figure 17 – Dates of all previous measurements for the corresponding patient


Figure 18 – Plaque report

You can save all the parameters by pressing the “Save Results” button xlsx files, as well the grayscale and RGB image by pressing the “Save Gray Plaque” and “Save RGB Plaque”, respectively.
 You can find details about the content of the files, on Appendix D.

Manual

This option gives you the opportunity to read this operation manual of the system.

Settings

This option gives you the opportunity to see your personal details and to change your password.

View and Edit personal details

By pressing the “View Profile” button, you can the profile page (Figure 19).

Figure 19 – Doctor’s information
On the top right of the page, there is an icon which when it is pressed, you are able to edit your personal information (Figure 20). Then you can press the save button (second button on the top left of the page) and a message confirms that you have made your changes correctly (Figure 21). In case of failure to fill in a mandatory field, the appropriate warning message will be displayed (Figure 3).


Figure 20 – Edit doctor’s information

Figure 21 – Changes on doctor's information has been changed

Change of Password


Here, you have the possibility to change your password. You have to fill the “New Password” and “Confirm” fields with your new password.

Figure 22 – Change Password

A screenshot of a web application interface for setting a password. The background is a light blue gradient with a stethoscope image on the right. At the top, a dark blue header bar contains a back arrow icon and the text "Set your password". In the center, there is a white rectangular form titled "Set password". Inside this form, there are two text input fields: "New Password" and "Confirm". Below these fields is a blue "Submit" button. A thin blue border highlights the "New Password" and "Confirm" fields.

The password must be the same in both fields to make the change. In case of the two fields are different, a warning message appears on the screen (Figure 23). Once the password is changed, you are moved to the main menu (Figure 9).

Figure 23 – Password must be the same in both fields

A screenshot of a warning message dialog box. It has a white background with a thin grey border. The text "BOTH FIELDS MUST BE THE SAME." is centered in a dark grey font. At the bottom right, there is a grey button with the text "OK" in a dark grey font.

Exit

This option disconnects you from your account, and returns back to the Login page (Figure 4).

Operations

This page includes all the main operations you can do (Figure 24).

Figure 24 – Main Operations


Image Histogram Normalization

Ultrasound (U/S) images can be normalized through this module. Image normalization is performed in order to be able to have reproducibility of texture features and plaque types (Elatrozy et al, 1998; Nicolaides et al, 2005). Blood and adventitia are used as two reference points and the image is linearly scaled so that zero is the gray scale value for blood and 190 for adventitia. In order to have images that can be normalized through the program in a meaningful way ultrasonographers should follow some basic but simple steps in image acquisition as described in Appendix A. First load an image that needs to be normalized (Figure 25).



Figure 25 – Image normalization: step 1.

User can press Zoom in/Zoom out if required. Then perform the following steps:

- 1) Press Define Low.
- 2) Select two points on the image to define the blood area. Try avoiding areas with “noise”.
- 3) Press Define High.
- 4) Select two points on the image to define the adventitia area (middle 2/4^{ths} of the brightest part).
- 5) Press Start Normalization

It is essential to use the zoom facility prior to step 3 in order to be able to select the middle 2/4^{ths} of the brightest section of adventitia adjacent to the plaque (Figure 26).

When normalization is finished, the normalized image is shown on the right site of the screen (Figure 27).


Figure 27 – Normalized image is shown to

Normalized image

The normalized image can be saved on the server by pressing the save file button. The image will be saved as (image_filename)_**norm**.(type_of_original_image) . When the file is saved the check box next to the save file button will appear as checked.

Image Scale Definition / Measurements

Several measurements such as intima-media thickness (IMT), plaque thickness, adventitia to adventitia distance or plaque area can be made through this module. This module can be used as a stand-alone option without using other modules. The main panel of this module is shown in Figure 10. On the left side of the panel the user is able to select the drive on which the images have been saved. In the list box near the drive selection the user is able to select the black and white image that will be opened. On the right side of the panel the user can select the corresponding color image that will be displayed next to the black and white image in order to help in the identification of the plaque outline. This is useful for plaques with black areas.


Figure 28– Image scale/ Measurements screen

If a colored image is chosen to upload at the Grayscale canvas, the image automatically transformed to grayscale.

This screen has two groups of buttons; the first group is the group of zoom and define image scale (Figure 29) and the second the group of measurements (Figure 32).


Figure 29– First group of buttons

When the Def.Im.Scale button is pressed the cursor becomes a cross. The user has to select two points on the image scale (Figure 30) and press enter. An input dialog asking about the distance of the two points in millimetres will appear for final scale definition (Figure 31). Then, the user has to press Submit scale.


Figure 30– Scale – Selection of points


Figure 31– Final Scale Definition

The next group of buttons is the measurements group (Figure 32). The image scale has to be defined before these buttons can be used

Figure 32 – U/S image Measurements

The first measurement is the IMT. After pressing this button the cursor becomes a cross and the user can select any number of pairs of points on the image and Submit IMT (Figure 33). The mean, max and min IMT is automatically calculated and displayed. If user selects an odd number of points, then the last number is ignored.


Figure 33 - Intima media Thickness

The second measurement of this group is the plaque thickness. When pressing this button the cursor becomes a cross and the user has to select a pair of points on the image and press Submit plaque (Figure 34).


Figure 34 – Plaque Thickness

The third measurement is the same as the second one. It is used in vessels that have two plaque components.

The fourth measurement is the vessel width (adventitia to adventitia). When pressing this button, the cursor becomes a cross and the user has to select as many pairs of points as he/she wants on the image and press Submit Adventitia (Figure 35). The mean vessel width is automatically calculated and displayed.


Figure 35– Vessel width

Finally, the area 1, 2 buttons are for measuring two different areas of the image (Figure 36). This is a free measurement tool that gives the selected area in mm².



Figure 36– Area

When finished with all actions the screen will appear.

Measurements are automatically saved in the database.

Change image visual resolution

Through this option the user is able to standardize whole or images of a previously cropped plaque to a fixed number of pixels per mm. Figure 37 shows the main screen of this option.

following.

1. Press the "Choose File" button.
2. Write "no" into the prompt question for uploading plaque or not.
3. Choose the image and click the "Upload" button.
4. Press the "Def.Im.Scale" button (see def. image scale / measurements section) to define the scale of the image. The cursor becomes a cross.
5. Select two points on the image scale and press enter.
6. Give the real distance of the two points in millimeters into the input dialog that is appeared.
7. Press "Change res 20 pix/mm button".

Note 1: The selected points must be on the same vertical line.

Note 2: The Zoom In/ Zoom Out action is very useful for the exact selection of the two points.

Note 3: During the step 5, the user is able to delete a point that was selected by mistake by pressing the 'd' button from the keyboard.

The resolution of the selected image will be changed to 20 pixels per mm according to the bicubic interpolation method and the new image will be saved as file_name_20res.(type_of_original_image).


Figure 38 – Input dialog to upload plaque

To change the resolution of a previously cropped plaque, (see image crop section) please do the following:

1. Press the “Choose File” button.
2. Write “yes” into the prompt question for uploading plaque or not.
3. Choose the image and click the “Upload” button.
4. Choose the new resolution and press the appropriate button “Change res 10pixels/mm, 15pixels/mm or 20pixels/mm” . The plaque file will be changed to the new resolution, using the bicubic interpolation method.

The

new file will be saved as

file_name_plq_num_10res.(type_of_original_image),

file_name_plq_num_15res.(type_of_original_image) or

file_name_plq_num_20res.(type_of_original_image) respectively.

Note: when cropping the plaque one should define the scale of the image. If this is not done then the resolution of the plaque image cannot be changed (see Image Crop section).


Figure 39 – Normalization of a plaque image

The images with the new resolution, are automatically saved into a folder named with the ID of the patient, on server.

In order to save the results to the database, the user must select the third icon at the top left of the page.

Image Crop

This module provides a means of selecting and cropping a plaque from an U/S carotid

image. Figure 40 shows the main panel of this module. On the left side the user is able to select the grayscale image that will be opened. On the right side of the panel the user may select the corresponding color image that can be displayed in order to help in the identification of the plaque outline (useful for plaques with black areas).

In order to change the contrast of the displayed image, the user can choose log image


Figure 40 – Image Crop main panel

The buttons (on the Select Action part and those for uploading the RGB image) cannot be selected until the grayscale image is selected.

For the image region crop, please do the following:

1. Upload the grayscale Image in order to let the button on the Select Action section to be enabled.
2. Upload the RGB image (not necessary).
3. Press the “Def.Im.Scale” button to define the scale of the image. The cursor becomes a cross →
4. Select two points on the image scale and press enter.
5. Give the real distance of the two points in millimeters into the input dialog that is appeared.
6. Press the “Image Region” button to select the region of interest. The cursor becomes a cross.
7. Select different points on the image in order to select the polygon outlining the plaque.
8. Press enter and the region will be displayed at the bottom of the panel.

Note 1: When an image, which is not gray scale such as RGB, or intensity image is loaded on the left screen, the program transforms the image automatically into Gray Scale.

Note 2: The user is able to change the contrast of the grayscale image in order to better distinguish the plaque area, by pressing the “Log Image” button.

Note 3: The Zoom In/ Zoom Out action is very useful for the exact selection of the points.

Note 4: When selecting points, the user is able to delete a point that was selected by mistake by pressing the ‘d’ button from the keyboard.

Note 5: If the “Image Region” button is selected without having defined the scale of the image, a warning message will be appeared (Figure 41).


Figure 41 – Warning message for Note 5.

Selection
of a polygon
on GrayScale
image


Figure 42 – Selections of polygon


Figure 43 – Cropped plaque appears

Explanation of step 7:

The user is able to select two plaque components (if necessary) from one image, as can be seen in the following figures (Figure 44 and Figure 45). In order to select the two plaques without having unnecessary pixels connecting them the user has to connect the two plaques by a thin line (Figure 45).


Figure 44 – Selection of plaques


Figure 45 – Plaques connected by a thin line

When the step 7 is done, the selected region is displayed on the “Image Region” canvas (Figure 46).


Figure 46 – Selected plaque

All the images among with the image at 20 pixels/mm which is automatically created, are saved into a folder named with the ID of the patient, on server. In order to save the results to the database, the user must select the third icon at the top left of the page.

Plaque texture parameters tool

The fifth module of the program is the plaque texture parameters tool. This function / tool can be selected either from the initial screen - features extraction (Figure 24) or from the image crop tool screen – feature extraction (Figure 47).


Figure 47 – Feature extraction

Firstly, the “Choose File” button must be selected and a popup window will be appeared (Figure 48). Choose the name of a plaque, and load this plaque. The system calculates the parameters as shown in Figure 49.


Figure 48 – Choose name of a plaque


Figure 49 – Parameters of plaque

Image contour coloring

In order to support plaque viewing, a coloring option has been added to this module.

The pixels of the cropped plaque are displayed in ten colors according to the following coloring scheme (Figure 50).

Gray levels	0-25: Black
	26-50: Blue
	51-75: Green
	76-100: Yellow
	101-125: Orange
	126-255: Red


Figure 50 – Colored plaque

Image texture features:

The texture features available make this module a powerful research tool. They consist of histogram measures, first order gray level parameters and Spatial Gray Level Dependence Matrices (SGLDM).

For definitions and brief descriptions of these texture features please see Appendix B.

Additional features option:

Through this option the user is able to record some extra characteristics of the plaque

such as the presence of a dark area adjacent to the lumen, patient status (Asymptomatic, AF, TIA, or Stroke), the presence of discrete white (red) areas in the plaque, type of plaque according to the Geroulakos classification and percentage of stenosis. The type of plaque is calculated automatically by the program, which quotes both plaque type 4 and 5 as 4. However, the plaque type can be changed manually.

Furthermore, predicted stroke risk rates, based on the results from the ACSRS study, can be calculated through this module. User has to press the ACSRS eq.risk button, answer the question about history of contralateral TIAs and stroke; and the predicted stroke risk rates will be calculated and displayed (Figure 52).

For details on how to calculate the stroke risk rates, please see Appendix C.


A screenshot of a software input dialog box. The dialog has a light gray border and a white background. At the top, the text "Previous history of Stroke and TIA? (yes/no)" is displayed in a small, black, sans-serif font. Below this text is a single-line text input field with a thin blue border. At the bottom of the dialog, there are two buttons: "OK" on the left and "Cancel" on the right. Both buttons have a light gray background and a thin blue border. The "OK" button is slightly larger and more prominent than the "Cancel" button.

Figure 51 - Input dialog for previous history of Stroke or TIA


Figure 52 – Predicted stroke risk rates based on the ACSRS results

Patient's Personal information

As can be seen from the pictures of the functions, once the doctor chooses or records a patient, then the patient's name appears at the top of the page on the right along with his / her identity (Figure 53). Next to the data, a button is displayed which, when selected, enables the physician to update and / or process the patient's data (Figure 5.49).

The physician can process all of the data, except for the patient's identity, for security issues. When a change is made, the save button is selected (second button on the top left), and the appropriate confirmation message appears

Figure 53 – Patients ID

Appendix A

Instructions to Ultrasonographers

Recommended settings of duplex scanner for plaque image analysis (El-Atrozy et al, 1998)

For carotid and femoral bifurcation scanning, a high-frequency linear array transducer(ideally 7 10 or 7-12 MHz) should be used. Several technical ultrasound settings need to be observed to ensure optimum image quality for measurements of IMT, plaque thickness and analysis of plaques / images.

1. Maximum dynamic range should be used. The dynamic range refers to the way in which gray scale information is compressed. A high dynamic range ensures a greater display of gray scale values which is important for ultrasonic tissue characterization.
2. Persistence should be set on low and frame rate set on high, the latter ensuring good temporal resolution
3. The Time Gain Compensation curve (TGC) should be sloping through the tissues but should be positioned vertically through the lumen of the vessel, because there is little attenuation of the beam as it passes through blood. This will ensure that the adventitia of the anterior wall has the same brightness as the adventitia of the posterior wall.
4. The overall gain should be adjusted to give optimum image quality. It should be turned up until there is noise in the lumen of the artery and then reduced to minimise but not abolish the noise. Areas of blood without noise should be available.
5. A linear post-processing curve should also be used.
6. Where possible the ultrasound beam should be at 90° to the arterial wall.

Measurements

1. IMT of the Common Carotid Artery (mm)

This is a measurement of the Intima-Media thickness (IMT) on the far wall of the common carotid artery 1-2 cm proximal to the bifurcation. It is well documented that measurements taken from the near wall are less accurate and therefore less reproducible than those taken from the far wall. This is because the adventitia is more echogenic than the media and bright echoes produced by the adventitia of the near wall can “spill” into the adjacent media . Thus, echoes from the media are lost. This effect is far less apparent on the far wall where the media and media-adventitia interface are closer to the probe than the adventitia.

2. Presence of Plaque

Plaque is defined as a localized thickening involving the intima and media in the bulb, internal carotid, external carotid or common femoral arteries which is greater than 1.2 mm.

APPENDIX B

Image texture features (Description of algorithms used)

Several parameters are available in order to provide the user a wide choice. Based on the gray level values the parameters are divided into several groups as shown in the plaque parameters tool.

First order gray level parameters

In this category [[3]] the parameters are derived from the gray level histogram. They describe the gray level distribution without considering spatial independence; as a result they can only describe echogeneity of texture and the overall variation characteristics within the Region Of Interest (ROI). The parameters used from this category are:

The Mean gray level (MEA):

$$MEA = \frac{1}{N} \sum_i \sum_j g(i,j)$$

this parameter describes the pixels mean gray level in a ROI.

Variance (VAR):

$$VAR = \frac{1}{N} \sum_i \sum_j (g(i,j) - MEA)^2$$

This parameter describes the histogram width or an estimation whether the differences of the pixels gray level in the ROI differs from the mean gray level.

Kurtosis (KUR):

$$KUR = \frac{1}{N} \left[\frac{\sum_i \sum_j (g(i,j) - \mu)^4}{\sigma^4} \right]$$

Median (MED): Is the median gray level of image gray levels

Skewness (SKE):

$$SKE = \frac{\left[\sum_i \sum_j (g(i,j) - \mu)^3 \right]}{N}$$

Energy (ENER) :

$$ENER = \sum_i p(i)^2$$

Entropy (ENT) :

$$ENT = -\sum_i p(i) \log(p(i))$$

where $g(i,j)$ is the gray level of the pixel (i,j) , N is the total number of pixels in the specified ROI, μ & σ are the mean value and standard deviation of total number of pixels, p is the probability distribution table.

The Spatial Gray Level Dependence Matrices (SGLDM)

The SGLDM [[8], [10]] algorithm is based on the assumption that texture properties of an image are contained in the overall or “average” spatial relationship between the gray levels in the image. SGLDM is based on the estimation of the second order conditional probability density $f(i,j;d,\theta)$. Each value $f(i,j;d,\theta)$ represents the probability that two different resolution cells which are in the direction specified by an angle θ and have distance d , will have gray level values i and j respectively.

Assume that we have to analyse a rectangular image of size $N_x \times N_y$ (where N_x are the resolution cells in the x direction and N_y in the y direction). Each cell's gray level value is quantitised to N_g levels. Let $L_x = \{1, 2, \dots, N_x\}$ be the horizontal spatial domain, $L_y = \{1, 2, \dots, N_y\}$ be the vertical spatial domain, $G = \{1, 2, \dots, N_g\}$ be the quantitised gray level values of the image. $L_x \times L_y$ is the set of resolution cells of the image. The image can be represented as a function $I(x,y)$, which has quantitised values from the set $1, \dots, N_g$ for each different pair of coordinates (x,y) , $I: L_x \times L_y \rightarrow G$.

The estimation of $f(i,j;d,\theta)$ is evaluated for angles quantitised to 45° . Normalised frequencies are defined by:

$$\begin{aligned} P(i, j, d, 0^\circ) &= \# \{ ((k,l), (m,n)) \in (L_y \times L_x) \\ &\quad \times (L_y \times L_x) : k - m = 0, |l - n| = d, \\ &\quad I(k,l) = i, I(m,n) = j \} \\ P(i, j, d, 45^\circ) &= \# \{ ((k,l), (m,n)) \in (L_y \times L_x) \\ &\quad \times (L_y \times L_x) : (k - m = d, l - n = -d) \\ &\quad \text{or } (k - m = -d, l - n = d), \\ &\quad I(k,l) = i, I(m,n) = j \} \\ P(i, j, d, 90^\circ) &= \# \{ ((k,l), (m,n)) \in (L_y \times L_x) \\ &\quad \times (L_y \times L_x) : |k - m| = d, l - n = 0, \\ &\quad I(k,l) = i, I(m,n) = j \} \\ P(i, j, d, 135^\circ) &= \# \{ ((k,l), (m,n)) \in (L_y \times L_x) \\ &\quad \times (L_y \times L_x) : (k - m = d, l - n = d) \\ &\quad \text{or } (k - m = -d, l - n = -d), \end{aligned}$$

$$l(k,l) = i, l(m,n)=j \}$$

Where # denotes the number of elements in the set.

If texture is coarse and d is small, compared to the texture elements, the pairs of points with distance d will have similar gray level values so the points on the main diagonal of the matrices $P(i,j;d,\theta)$ will have high values. If texture is smooth the values will be more spread out.

The characterisation is done using features which measure the scatter around the main diagonal of the matrices. Harallick et al. (1973) has proposed 14 measures to extract useful texture information from an image. Some of the measures relate to specific textural characteristics of the image, such as homogeneity, contrast, the presence of organised structure within the image, complexity and nature of gray tone transitions.

Despite the fact that several characteristics contain information about certain textural features, it is difficult to identify which characteristic is identified by each measure.

The program calculates the following features:

Angular Second Moment (ASM):

$$f_1 = \sum_i \sum_j \{p(i,j)\}^2$$

Contrast (CON):

$$f_2 = \sum_{i=1}^{N_g-1} n^2 \left\{ \sum_{i=1}^{N_g} \sum_{j=1}^{N_g} p(i,j) \right\}_{|i-j|=k}$$

Correlation (COR):

$$f_3 = \frac{\sum_i \sum_j (ij) p(i,j) - \mu_x \mu_y}{\sigma_x \sigma_y}$$

Sum of squares: Variance (VAR):

$$f_4 = \sum_i \sum_j (i-j)^2 p(i,j)$$

Inverse Difference Moment – Homogeneity (HOM):

$$f_5 = \sum_i \sum_j \frac{1}{1+(i-j)^2} p(i,j)$$

Sum Average (SAV):

$$f_6 = \sum_{i=2}^{2N_g} i p_{x+y}(i)$$

Sum Variance (SVA):

$$f_7 = \sum_{i=2}^{2N_g} (i - f_8)^2 p_{x+y}(i)$$

Sum Entropy (SEN):

$$f_8 = - \sum_{i=2}^{2N_g} p_{x+y}(i) \log\{p_{x+y}(i)\}$$

Entropy (ENT):

$$f_9 = - \sum_i \sum_j p(i, j) \log(p(i, j))$$

Difference Variance (DVA):

$$f_{10} = \text{variance of } p_{x-y}$$

Difference Entropy (DEN):

$$f_{11} = - \sum_{i=0}^{N_g-1} p_{x-y}(i) \log\{p_{x-y}(i)\}$$

Information Measures of Correlation:

(InM1):

$$f_{12} = \frac{HXY - HXY1}{\max\{HX, HY\}}$$

(InM2):

$$f_{13} = (1 - \exp[-2.0(HXY2 - HXY)])^{1/2}$$

$$HXY = - \sum_i \sum_j p(i, j) \log(p(i, j))$$

Where HX and HY are the entropies of p_x, p_y and

$$HXY1 = - \sum_i \sum_j p(i, j) \log\{p_x(i) p_y(j)\}$$

$$HXY2 = - \sum_i \sum_j p_x(i) p_y(j) \log\{p_x(i) p_y(j)\}$$

where N_g is number of gray levels, μ_x & σ_x are the mean and standard deviation of the row sums of the matrix $p(i, j)$, and μ_y & σ_y are the corresponding statistics of the column sums with p_{x+y} and p_{x-y} given by:

$$p_{x-y}(k) = \sum_{i=1}^{N_g} \sum_{j=1}^{N_g} p(i, j), \quad k=0,1,\dots,N_g-1$$

$$|i-j|=k$$

$$p_{x+y}(k) = \sum_{i=1}^{N_g} \sum_{j=1}^{N_g} p(i, j), \quad k=2,3,\dots,2N_g$$

$$i+j=k$$

Each measure is evaluated for $d=10^1$ and $\theta=0^\circ, 45^\circ, 90^\circ, 135^\circ$; features are obtained from the sample mean of the estimated values.

Appendix C

Calculation of predicted 5 year survival based on the ACSRS study results

The hazard ratios for covariates presented in table 1 can be combined with the baseline survival function to predict 5-year survival probabilities for a patient with a given set of covariates.

For an individual with baseline covariate values (close to mean covariates) as listed in table 2 the stroke free survival at 5 years, $S_0(5y)$ is estimated as 0.972.

Table 1

Proportional hazards models with ipsilateral CORI events as the dependent variable.

(GSM = Gray Scale Median; DWA = Discrete White Areas; HR = Hazard Ratio).

Clinical factors with plaque features. Ipsilateral hemispheric stroke as the dependent variable. 5-year baseline hazard estimated as .972; Harrell's C = .80; Pseudo R² = .61

Variable	β	HR	95% CI	P
Stenosis (% ECST)	0.026	1.03	1.00 to 1.05	0.024
Log(GSM + 40)	-2.672	0.07	0.02 to 0.20	<0.0001
Plaque area ^{1/3} (mm ³)	0.629	1.88	1.28 to 2.75	<0.0001
DWA (Present vs. absent)	0.429	1.54	0.81 to 2.92	0.18
History of contralateral TIAs and/or stroke (Yes vs. no)	0.973	2.65	1.54 to 4.54	<0.0001

Table 2. Baseline covariates and transformed values in an individual

Covariate	Baseline value	Corresponding transformations
Stenosis, %	80	None
GSM	30	$\ln(\text{GSM}+40) = 4.248$
Plaque area, mm ²	40mm ²	$(\text{Plaque area})^{1/3} = 3.420$
DWA	Absent	0
History of contralateral TIAs and/or stroke	Absent	0

Predicted 5-year stroke free survival for a patient with different values of the covariates in table 2 can be calculated as follows.

1. Transform continuous covariates using the formulae given in table 2, third column.
2. Using these (possibly transformed) values, calculate the differences, x , between the value you are interested in and the values used in $S_0(5y)$ (Table 3, second column).
3. Multiply each of these differences by the corresponding log-hazard ratio, β , from table 1 to obtain βx as shown in table 3 column 4.
4. Sum the values obtained in step 3 above to denote this as βx .

5. Calculate $\exp(\beta x)$
6. Compute predicted survival probability as a percentage for a patient using $100 \times 0.972^{\exp(\beta x)}$.

Following this calculation for a patient with 90% stenosis, GSM = 10, plaque area = 80 mm², DWA present and history of contralateral TIAs and/or stroke absent:

1. $\ln(\text{GSM}+40) = 3.912$. $(\text{Plaque area})^{1/3} = 4.309$

Table 3 Steps 2 and 3 in the calculation of 5 year predicted stroke free survival

Variable	Calculation	2. x	β	3. βx
Stenosis	$90 - 80$	10	0.026	0.2600
$\ln(\text{GSM}+40)$	$3.912 - 4.248$	-0.336	-2.672	0.8978
$(\text{Plaque area})^{1/3}$	$4.309 - 3.4200$	0.889	0.629	0.5592
DWA	$1 - 0$	1	0.429	0.4290
History of contralateral TIAs and/or stroke	$0 - 0$	0	0.973	0

4. Sum of $\beta x = \beta x = 0.2600 + 0.8978 + 0.5592 + 0.4290 + 0 = 2.1460$
5. $\exp(2.146) = 8.551$
6. Predicted 5-year stroke free survival $100 \times 0.972^{\exp(2.146)} = 78.4\%$

Appendix D

Description of database tables

There are eight tables in the database, which are explained below:

- **Administrators Data table (adminData):**

This table includes the details of system administrators, who can add a new physician.

Specifically, this table includes the following:

- id: This field is the primary key of the table and is automatically generated for each administrator in ascending order.
- email: The administrator email, varchar type. Obligatory field.
- password: This is the code of the administrator to be able to connect to the system, encrypted using the Secure Hash Algorithm 256.

- **Patients Data Table (patientData):**

This table includes the demographics of the patient so the doctor can identify him when he wants to process his own ultrasound images.

Specifically, this table includes the following:

- id: This field is the primary key of the table and is automatically generated for each patient in ascending order.
- patient_id: It is the patient's identity which is completed by the doctor. It is a varchar type. The field is mandatory.
- national_healthcare_patient_id: It is the number of the patient's international medical card. It is a varchar type. Non-mandatory field.
- given_name: The name of the patient, type varchar. Obligatory field.
- family_name: The patient's surname, varchar type. Obligatory field.
- date_of_birth: The date of birth of the patient, type date. Obligatory field.
- gender_coding: The patient's gender. M - Male, F - Female. Obligatory field.
- telephone_no: Patient phone number, int type. Obligatory field.
- email: Patient e-mail, varchar type. Non-mandatory field.

- address_id: The key that appears in the addressData table, and specifies the patient address details (Foreign Key).
-
- Doctors Data Table (doctorData):
 This table includes the demographics of the doctor, which the administrator states when he adds a new doctor. The doctor himself can then process it.
 In particular, it includes the following fields:
 - di: This field is the primary key of the table and is automatically generated for each doctor in ascending order.
 - doctor_id: It is the doctor's identity. It is a varchar type. The field is mandatory.
 - national_healthcare_doctor_id: It is the number of the doctor's international medical card. It is a varchar type. Non-mandatory field.
 - given_name: The name of the doctor, type varchar. Obligatory field.
 - family_name: The doctor's surname, varchar type. Obligatory field.
 - date_of_birth: The date of birth of the doctor, type date. Obligatory field.
 - gender_coding: The doctor's gender. M - Male, F - Female. Obligatory field.
 - telephone_no: Doctor phone number, int type. Obligatory field.
 - email: Doctor e-mail, varchar type. Mandatory field.
 - pass_word: It is the code that the physician declares to be able to connect to the system, encrypted using the Secure Hash Algorithm 256.
 - address_id: The key that appears in the addressData table, and specifies the doctor address details (Foreign Key).

- Physician – patient table (doctor_patient):

In this table there are the codes of doctors and their patients, so as to distinguish which doctor can access each patient's data and process his data.

In particular, it includes the following fields:

- di: The primary key of the table, automatically generated by the system, in ascending order.
- doctor_id: Foreign Key, corresponds to the primary key of the doctorData table, and states the doctor.
- patient_id: Foreign Key, corresponds to the primary key of the patientData table and indicates the patient.

- Addresses Data Table (addressData):

In this table can be found the details of the address of both the patient and the doctor.

Specifically, the address fields are:

- address_id: The primary key, which is automatically generated in ascending order, is what corresponds to the patient / physician with the addresses.
- street: Corresponds to street name. Field type varchar.
- street_no: Address number. Field type Int.
- province: Corresponds to the province or area of residence. Field type varchar.
- city: The city where the patient / doctor resides. Varchar field.
- country: The country of residence. Field type varchar.
- post_code: Postcodes. Varchar field.

- Table with the data of the ultrasound images (ultrasoundImages):

Table of data storage of each ultrasound image of the patient's carotid artery.

In particular, the following elements are included:


- ci: The primary key of the table, automatically generated by the system, in ascending order.
 - patient_id: the patient's mapping key of the ultrasonic image. It is matched to the primary key of the patientData table.
 - date: The date when this image was processed.
 - time: The time at which this image was processed.
 - gray_image: The ultrasonic image in grayscale form.
 - normalized_image: The normalized image created when performing the first operation.
 - rgb_image: The ultrasonic image in color.
 - res20_image: The ultrasonic image with a 20 pixels / mm resolution generated during the third and fourth operations.
 - scale: The scale of the image as stated by the doctor.
-
- Table with the IMT parameters (IMT_Parameters):
This table includes all the measurements related to IMT.
Specifically, this table includes:
 - ci: The primary key of the table, automatically generated by the system, in ascending order.
 - patient_id: the patient's mapping key of the ultrasonic image. It is matched to the primary key of the patientData table.
 - date: The date when this image was processed.
 - time: H The time at which this image was processed.
 - ppmm: The pixels per millimeter of the image scaling.
 - imtmax: The maximum number of IMT thickness.
 - imtmin: The minimum number of IMT thickness.
 - imtmean: The mean value of IMT thickness.
 - advmean: The mean value of vessel width.

- plq1: The thickness of plaque 1.
 - plq2: The thickness of plaque 2.
 - area1: The area of plaque 1.
 - area2: The area of plaque 2.
 - gray_image: The image in grayscale form.
 - rgb_image: The image in RGB form.
- Table with the elements and parameters of the plaque (plaqueParameters):
The last table of the database contains the data and measurements taken during the carotid plaque analysis.
The fields in this table are as follows:
 - ci: The primary key of the table, automatically generated by the system, in ascending order.
 - patient_id: the patient's mapping key of the ultrasonic image. It is matched to the primary key of the patientData table.
 - ultrasound_id: the plaque's mapping key of the ultrasonic image. It is matched to the primary key of the ultrasoundImages table.
 - date: The date when this image was processed.
 - time: H The time at which this image was processed.
 - plaque: The region of plaque on the ultrasonic image.
 - mask: Mask created by the system, which has the value 1 at the points corresponding to the foreground of the plaque image, and the value 0 at the points corresponding to the background of the plaque. It is a necessary field for calculating the features only from the foreground points of the plaque.
 - plaque_10res: Plaque image with new resolution 10pixels / mm created in operation 3.
 - plaque_15res: Plaque image with new resolution 15pixels / mm created in operation 3.

- plaque_20res: Plaque image with new resolution 20pixels / mm created in operation 3.
- kurtosis: The value of kurtosis calculated by the system.
- total_num_pixels: The total number of pixels that make up the foreground of the plaque image.
- mm2: The area of the plaque in square millimeters.
- below30: The pixel percentage of the plaque having a value less than 30.
- below50: The pixel percentage of plaque having a value greater than 30 but less than 50.
- mean: The average of the pixel values contained in the plaque.
- variance: The variance of the plaque.
- median: The median value from the pixel values of the plaque.
- skewness: The skewness value.
- energy: The energy value.
- entropy: The entropy value.
- haral_ASM: The value of the angular second moment by the Haralick method.
- haral_contrast: The value of the contrast by the Haralick method.
- haral_correlation: The value of the correlation by the Haralick method.
- haral_variance: The value of the variance by the Haralick method.
- haral_homogeneity: The value of the homogeneity by the Haralick method.
- haral_sum_average: The value of the sum average by the Haralick method.
- haral_sum_variance: The value of the sum variance by the Haralick method.
- haral_sum_entropy: The value of the sum entropy by the Haralick method.
- haral_entropy: The value of the entropy by the Haralick method.
- haral_diff_var: The value of the difference of variance by the Haralick method.

- haral_diff_entr: The value of the difference of entropy by the Haralick method.
- haral_inf_meas1: The value of the Information Measures of Correlation 1 by the Haralick method.
- haral_inf_meas2: The value of the Information Measures of Correlation 2 by the Haralick method.
- scaling: The scale of the image as stated by the doctor.
- darkArea: Indication of presence of a dark area next to the lumen. 'Y' - if any, 'N' - if it does not exist.
- status: Patient status indication ('A' - Asymptomatic, 'B' - AF, 'C' - TIA, 'D' - Stroke)
- dwa: Indication of presence of distinct white areas on the plaque. 'Y'- if they exist, 'N' - if they do not exist.
- type: Plaque type ('1', '2', '3', '4', '5').
- perc: The stenosis rate of the plaque.
- previous_history: Indication of patient history. 'Y' - if history exists, 'N' - if there is no history.
- year_5_stroke: The survival rate of the patient for 5 years without a stroke.
- aver_annual_stroke: The average annual rate of stroke.
- lnGSM: The value $\ln(\text{GSM} + 40)$ calculated by the system, where GSM is the Gray Scale Median.
- area13: The value $(\text{Plaque area})^{1/3}$ calculated by the system.
- bel10: The number of pixels in the plaque having a value of less than 10.
- bel20: The number of pixels in the plaque having a value of less than 20 and greater than 10.
- bel30: The number of pixels in the plaque having a value of less than 30 and greater than 20.
- bel40: The number of pixels in the plaque having a value of less than 40 and greater than 30.
- bel50: The number of pixels in the plaque having a value of less than 50 and greater than 40.

- col1: The number of plaque pixels having a gray level value of less than 25.
- col2: The number of plaque pixels having a gray level value of less than 50 and greater than 25.
- col3: The number of plaque pixels having a gray level value of less than 75 and greater than 50.
- col4: The number of plaque pixels having a gray level value of less than 100 and greater than 75.
- col5: The number of plaque pixels having a gray level value of less than 125 and greater than 100.
- col6: The number of plaque pixels having a gray level value of less than 150 and greater than 125.
- col7: The number of plaque pixels having a gray level value of less than 175 and greater than 150.
- col8: The number of plaque pixels having a gray level value of less than 200 and greater than 175.
- col9: The number of plaque pixels having a gray level value of less than 225 and greater than 200.
- col10: The number of plaque pixels having a gray level value of less than 255 and greater than 225.

Appendix E

Description of Files Content

The content of IMT Report is :

	A	B	C	D	E	F	G	H	I	J
1	Patient Name	Plaque Name	Pixels per mm	Maximum Intima Media Thickness	Minimum Intima Media Thickness	Mean Intima Media Thickness	Mean Adventitia to Adventitia	Plaque 1	Plaque 2	Area 1
2	Maria Georgiou , 101010	15rc2bwn.tif	16,1	0,559	0,248	0,408	5,245	2,095	2,874	290,093
3										

The content of Plaque Report is:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Patient Na	Plaque Na	pixels<10	10<pixels<20	20<pixels<30	30<pixels<40	40<pixels<Below 30	Below 50	Col1	Col2	Col3	Col4	Col5	Col6	Col7	
2	Andreas Sz	15rc2bwn	3681	3178	2687	2338	1627	50,78	21	8573	5116	3228	1219	494	170	0
3																

The fields are explained at the Appendix D, in the corresponding tables of database.

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Technical Manual

Contents

General Instructions for the installation of the program on server	1
Run the Program	2
Modify the Database	3
Add a new Doctor	3
Forgot Password	4
HTML Files	4
View's Functions	5
Database	11

General Instructions for the installation of the program on server

For the transformation of the system onto the server, the following must be done:

1. If the settings.py file has something else than those starting with "django.xxx" and the name of the app on INSTALLED_APPS, then the user must install it with the pip.
2. In wsgi.py, the user must write the following code:

```
import os, sys
sys.path.append(" ") #the path of the app
```

```
os.environ.setdefault("DJANGO_SETTINGS_MODULE", " ") #the name of the
folder where the settigs.py is located
```

```
from django.core.wsgi import get_wsgi_application
application = get_wsgi_application()
```

3. Then, he have to go to server on path: /etc/apache2/apache2.conf, and write the following code on line 175:

```
WSGIScriptAlias / ..... # the path of wsgi.py file
WSGIProxyPath ..... # the path of the app
```

```
Alias ..... #The path where to collect the static files
```

```
<Directory ..... >#The path where to collect the static files
Require all granted
</Directory> #The user gives access on this path to the apache
```


```
<Directory > # the path of the app
<Files wsgi.py>
Require all granted
</Files>
</Directory> #The user gives access on wsgi.py to the apache
```

4. Restart the apache service: `sudo service apache2 restart`
5. Run this command on terminal: `python3 manage.py runserver 173.249.10.105:8000`
6. If errors occur, then all the errors can be seen on `/var/log/apache2/error.log`.

Run the Program

In order to run the program, you must login into the server(Figure 1), with the following details:

- IP: 173.249.18.105
- Port Number: 22
- Username: ucy
- Password: U%c57!&kl

The server must be running during the usage of the program (Figure 2). For this, you must open a terminal and change the directory in thesis2018 using this command: “cd thesis2018/”. Then write the following command:

- python3 manage.py runserver 173.249.10.105:8000

```

ucy@vmi154301: ~/thesis2018
Using username "ucy".
ucy@173.249.18.105's password:
Welcome to Ubuntu 16.04.4 LTS (GNU/Linux 4.4.0-103-generic x86_64)

 * Documentation:  https://help.ubuntu.com
 * Management:    https://landscape.canonical.com
 * Support:       https://ubuntu.com/advantage

  _____
 |  _   _  |
 | | | | | |
 | |_| | | |
 |  _  | | |
 | | | | | |
 |_|_|_|_|_|

Welcome!

This server is hosted by Contabo. If you have any questions or need help,
please don't hesitate to contact us at support@contabo.com.

Last login: Sat Jun  2 19:14:56 2018 from 62.228.130.60
ucy@vmi154301:~$ cd thesis2018/
ucy@vmi154301:~/thesis2018$
ucy@vmi154301:~/thesis2018$ python3 manage.py runserver 173.249.18.105:8000
Performing system checks...

System check identified no issues (0 silenced).
June 02, 2018 - 17:16:20
Django version 2.0.5, using settings 'thesis.settings'
Starting development server at http://173.249.18.105:8000/
Quit the server with CONTROL-C.
```

Modify the Database

Models.py file:

This file contains the database tables of the system. If the user wants to make some changes in database, then he can make the changes into this file. Then, he must run the following command in terminal:

- python3 manage.py makemigrations
- python3 manage.py migrate

Add a new Doctor

The administrators of the program are the only ones who have the ability to add new doctors, for security issues.

For the moment, there is one administrator in the program, who can login with the following:

Email: epl400ADE2018@gmail.com

Password: 1234

URL: <http://173.249.18.105:8000/admin> login/


Forgot Password

If a doctor forgets his password, then an email is sent to him in his personal email address, with a new random password. This email is sent via the administrator of the program, namely via the mail: epl400ADE2018@gmail.com

To login the account of the administrator the password is: thesis2018.

HTML Files

The following diagram aims to show how the html files are connected in the system.

View's Functions

All the forms that is submitted and every request calls a function in the views.py. In this section you can find a brief description of every function:

login:

When requested, this function brings user to the login.html page. When the user submits his email and password, function saves the id into session 'doctor_id', in order to use it later. Function encodes the password with sha256 and compares the digest with the digest saved in the database. If the digests match, then user proceeds to menu.html.

admin_login:

When requested, this function brings admin to the admin_login.html page. When the admin submits his email and password, function encodes the password with sha256 and compares the digest with the digest saved in the database. If the digests match, admin proceeds to add_new_doctor.html.

edit_profile:

When requested, this function brings user to the edit_profile.html page. When user requests to change information, fields are enabled for change and function

gets from the page the content of all fields and updates database tables: doctorData.

Edit_patient_profile:

When requested, this function brings user to the edit_patient_profile.html page. When user requests to change information, fields are enabled for change and function gets from the page the content of all fields and updates database patientData.

Profile:

When requested, this function brings user to the my_profile.html page. Function loads from database all information from table doctorData and addressData

Id_generator:

This function creates a random sequence of 8 characters for system to use as new password, in case user forgets his.

Change_pass:

When requested, this function brings user to the change_pass.html page. When user submits new password, function checks if both fields are the same. If not, then function returns at same page with a message, warning user to use same password in both fields. Else, function encodes new password with sha256 and updates database table doctorData.

Forgot_pass:

When requested, this function brings user to the forgot_pass.html page. When user submits new password, function checks email exists in database. If not, the function returns at the same page. Else, function calls id_generator() and via the email account epl400ADE2018@gmail.com sends the product of id_generator() to user, and updates database table doctorData.

Add_new_doctor:

When requested, this function brings user to the add_new_doctor.html page. When user submits form, function adds new record in database table patientData.

Add_new_patient:

When requested, this function brings user to the add_new_patient.html page. Function, checks if patient id already exists in database. If not, then function adds new record in database table patientData. Else, if patient exists and was added to the system by the same user, then function does not allow user to add patient again. If patient exists, but was added by a different user, then allows user to use patient's existing information.

Operations:

When requested, this function brings user to the Operationion.html page. Function also sends patients full name and id, so they can be displayed at the top-right part of the page.

Menu:

When requested, this function brings user to the menu.html page. Function checks if doctor_id is null, meaning that no one signed in, and it returns to the login. If user selects "browse patient" from option "Patients" then function returns a list of all patients registered by the specific user. When user selects patient, then function saves patient's id, name, last name, table's primary key and plaque id into sessions to use it later. When user selects "Exit", session doctor_id becomes null.

Report:

When requested, this function brings user to the report.html page. This function creates a list with all dates of previous IMT measurements of a specific patient. When the date is chosen, function brings user at report.html page. If user chooses "Save Results", then an excel file is created and then downloaded to user's personal computer.

Imghist:

When requested, this function brings user to the imghist.html page. When user submits an image file, function saves image in folder '/static/' and shows image's histogram. When user submits low area, function gets the coordinates submitted and creates a sub image. After this, function calculates the median of low sub image. Some actions are done for submitting high area. With start_normalization, function uses previous low and high medians and applies image normalization to the original image. Finally, function saves normalized image and histogram and uploads it to imghist.html.

norm_res:

This function corresponds to the 3rd operation and interacts with the norm_res.html page.

Firstly, it checked if the doctor_id is empty, then no one has logged in to the system, so it returns to the login page. Otherwise, it does the following:

If the user wants to process a plaque image, then this function finds all the available plaques for the specific patient and sends them back to the html page. The user chooses the plaque he wants and accordingly to what number pixels per millimeter he wants to make the new image, the system creates that image.

If the user wants to process an ultrasonic image, then he must define the scale of the image, and using the Euclidean distance, the function calculates the pixels per millimeter of the image. It creates a new image with the new resolution (20pix/mm) and this image is displayed on the html page.

image_crop:

This function corresponds to the 4th operation, namely the Image Crop and interacts with the image_crop.html file. Firstly, it checked if the doctor_id is empty, then no one has logged in to the system, so it returns to the login page. Otherwise, it does the following:

When user wants to upload the grayscale or the RGB image, via image_crop.html, this function helps the uploading. It finds the location of images, and it send them back to the html file. In addition, it can log the grayscale image if the user wants to see the contrast more enhanced.

Also, in this function the scaling of the image is calculated in the same way as in the previous functions. Automatically, it generates the image with 20pix/mm in order to save it to the database.

If the user wants to save the data to the database, then this function calls another function (crop_plaque) which is responsible to crop the region of interest of the whole image. Then all the information is saved to the database.

crop_plaque:

This method is responsible for finding the appropriate polygon on the original grayscale image that corresponds the region of plaque, according to the points that user noted on the html page. It crops the plaque and creates a mask of the plaque. In the mask the pixels have the value of 1 or 0. The zero value on the

pixels that are inside the plaque region, and the one value on the pixels out of bounds of the plaque.

@points: The points (x,y) that user noted on the html file

@imgName: The name of the original grayscale image

@ppmm_plaque: The pixels per millimeter of the image

@top: The y-coordinate of the upper point of image

@left: The x-coordinate of the more left point of image

plaque_texture:

This function corresponds to the 5th module, and interacts with the plaque_texture.html file.

Firstly, it checked if the doctor_id is empty, then no one has logged in to the system, so it returns to the login page. Otherwise, it does the following:

It finds all the available cropped plaques for this patient and send them to the html file, in order to choose the doctor which plaque wants to process. Then, it calculates all the parameters, histogram measurements, first order statistics, SGLDM statistics and so on, for the specific plaque and saves them to the database.

When the physician wants to see the risk rates, then this function calculates the two risk rates and send them back to the html file.

find_average:

This function finds the average value of an array of values.

@x: The array of values

PolygonArea2:

This function calculates the area of a polygon region, specifically of the plaque region.

@x: The x-coordinates

@y: The y-coordinates

@ppmm: The pixels per millimeter of the image

createMask:

It creates a mask for the plaque, if there is no matching mask.

@maskname: The name for the new mask

@imgname: The name of the image

@direct: The path to save this mask

save_plaque_features:

This function takes the input parameters of the plaque, and calls another function to save it to the database.

saveResults:

This function save the input parameters of the plaque to the database.

@darkArea: Existence of dark area next to the lumen

@status: Patient status

@dwa: Existence of discrete white areas

@type: Plaque type

@perc: The percentage of stenosis

Plaque_report:

This function finds the parameters of plaque, according to the specific patient, date and time that were selected from the physician and sends them back to the Plaque_report.html page. Then, it creates an excel file with all these parameters and saves it on the local computer of the physician. Also, it saves the grayscale plaque and the rgb plaque.

Measurements:

When requested, this function brings user to the measurements.html page. When user submits an image file, function saves image in folder and updates database table utltradounsImages with the path. Some procedure is done if user uploads a colored image.

Def_im_scale_measurements:

When requested, this function brings user to the measurements.html page. Function takes the coordinates of two points and calculates the Euclidian distance between them. Also, function takes the number of millimeters and calculates pixels per millimeters(ppmm). Finally, function adds a new record at table IMT_Parameters, with the ppmm and the paths saved in function measurements. In session "IMT_id" the primary key of the record created id saved, in order to add the measurements of IMT that will follow.

IMT:

When requested, this function brings user to the measurements.html page. Function takes two lists, one with all x coordinates and one with all y coordinates. If the lists contain odd number of numbers, then the last number is ignored. Function calculates the Euclidian distance between them and with the ppmm found at Def_im_scale_measurements it finds the maximum distance, the minimum and the average. Finally, function adds the max, min and average value into table IMT_Parameters, at record saved in session "IMT_id".

Adventitia:

When requested, this function brings user to the measurements.html page. Function takes two lists, one with all x coordinates and one with all y coordinates. If the lists contain odd number of numbers, then the last number is ignored. Function calculates the Euclidian distance between them and with the ppmm found at Def_im_scale_measurements it finds the average distance. Finally, function adds average value of vessel width into table IMT_Parameters, at record saved in session "IMT_id".

Plaque_th:

When requested, this function brings user to the measurements.html page. Function takes the coordinates of two points and calculates the Euclidian distance between them and with the ppmm found at Def_im_scale_measurements it finds the real distance in mm, which represents the thickness of the plaque. Finally, function adds average value into table IMT_Parameters, at record saved in session "IMT_id".

Plaque_th2:

Same as plaque_th, when requested, this function brings user to the measurements.html page. Function takes the coordinates of two points and calculates the Euclidian distance between them and with the ppm found at Def_im_scale_measurements it finds the real distance in mm, which represents the thickness of the plaque. Finally, function adds average value into table IMT_Parameters, at record saved in session "IMT_id".

Area1:

When requested, this function brings user to the measurements.html page. Function takes two lists, one with all x coordinates and one with all y coordinates and with the function PolygonArea it calculates the area. The result is divided with the ppm found at Def_im_scale_measurements it finds the real area in mm. Finally, function adds area into table IMT_Parameters, at record saved in session "IMT_id".

Area2:

Same as area1, when requested, this function brings user to the measurements.html page. Function takes two lists, one with all x coordinates and one with all y coordinates and with the function PolygonArea it calculates the area. The result is divided with the ppm found at Def_im_scale_measurements it finds the real area in mm. Finally, function adds area into table IMT_Parameters, at record saved in session "IMT_id".

PolygonArea1:

This function calculates the area of a polygon region.

Database

The following diagram shows all the tables in the database and how they are connected.

The name and password are:

Username: root

Password: root